illumina®

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
※1 Lindeman NI, Cagle PT, Aisner DL, et al. Updated Molecular Testing Guideline for the Selection of Lung Cancer Patients for Treatment With Targeted Tyrosine Kinase Inhibitors: Guideline From the College of American Pathologists, the International Association for the Study of Lung Cancer, and the Association for Molecular Pathology. *J Mol Diagn.* 2018 Mar;20(2):129-159.

※2 Belli C, Penault-Llorca F, Ladanyi M, et al. ESMO recommendations on the standard methods to detect RET fusions and mutations in daily practice and clinical research. *Ann Oncol.* 2021 Mar;32(3):337-350.

※3 Chakravarty D, Johnson A, Sklar J, et al. Somatic Genomic Testing in Patients with Metastatic or Advanced Cancer: ASCO Provisional Clinical Opinion. *Journal of Clinical Oncology* 2022 40:11, 1231-1258

※4 Major US guidelines for NSCLC and CRC

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Genomic Findings with
**Evidence of
Clinical Significance**

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Genomic Findings with
**Potential
Clinical Significance**

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illumina | TruSight™ Oncology Comprehensive (JP) FOR IN VITRO DIAGNOSTIC USE Report Date 2023-02-07

Sample ID DNA5	Run QC	✓ PASS	Run ID 230206_NDX550320_0012_AHNYGGBDXX
Tumor Type Other solid tumor	RNA External Control & NTC	N/A	Analysis Date 2023-02-07
Sex Unknown	RNA Library QC	N/A	Knowledge Base Version 8.1.0.0
	DNA External Control & NTC	✓ PASS	Knowledge Base Published Date 2022-11-01
	DNA Library QC	✓ PASS	Module Version 2.5.1.55
	L DNA MSI QC	✓ PASS	Claims Package Version 1.1.0.0
	L DNA Small Variant & TMB QC	✓ PASS	
	L DNA Copy Number Variant QC	✓ PASS	

Alterations and Biomarkers Identified

The genomic findings reported below, for variants or biomarkers identified in this sample, are intended to provide tumor profiling information in accordance with professional guidelines.

Genomic Findings with Evidence of Clinical Significance *

Detected Variants	Details
BRAF p.(Val600Glu)	Type: SNV VAF: 10.86% Consequence: Missense Variant Nucleotide Change: NM_004333.6:c.1799T>A Genomic Position: chr7:140453136 Reference Allele: A Alternate Allele: T

Genomic Findings with Potential Clinical Significance *

TMB: 6.3 Mut/Mb		MSI: MS-Stable	
Detected Variants	Details		
AKT1 p.(Glu17Lys)	Type: SNV VAF: 7.8% Consequence: Missense Variant, Splice Region Variant Nucleotide Change: NM_001014432.2:c.49G>A; NM_001382430.1:c.49G>A Genomic Position: chr14:105246551 Reference Allele: C Alternate Allele: T		
APC p.(Arg1450Ter)	Type: SNV VAF: 10.08% Consequence: Stop Gained Nucleotide Change: NM_000038.6:c.4348C>T Genomic Position: chr5:112175639 Reference Allele: C Alternate Allele: T		
APC p.(Thr1556AsnfsTer3)	Type: Insertion VAF: 7.98% Consequence: Frameshift Variant Nucleotide Change: NM_000038.6:c.4666dup Genomic Position: chr5:112175951 Reference Allele: G Alternate Allele: GA		
ATM p.(Cys353SerfsTer5)	Type: Deletion VAF: 6.44% Consequence: Frameshift Variant Nucleotide Change: NM_000051.4:c.1058_1059del Genomic Position: chr11:108117845 Reference Allele: CTG Alternate Allele: C		
CTNNB1 p.(Thr41Ala)			

*Additional information in Informatics Details section

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EGFR p.(Leu858Arg)
EGFR p.(Thr790Met)

Type: SNV

VAF: 9.98% | Consequence: Missense Variant | Nucleotide Change: NM_005228.3:c.2573T>G | Genomic Position: chr7:55259515 | Reference Allele: T | Alternate Allele: G

Type: SNV

VAF: 3.88% | Consequence: Missense Variant | Nucleotide Change: NM_005228.3:c.2369C>G | Genomic Position: chr7:55249071 | Reference Allele: C | Alternate Allele: T

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